

HIV-1₄₅₁ gp120

ABL, INC.



Catalog Number:	5424
Quantity:	115 µg
Source:	Purified by affinity chromatography and HPLC, from inactivated, concentrated, conditioned medium of HUT78 (clone 6D5) cells chronically-infected with HIV-1 ₄₅₁ isolate.
Inactivation Procedure:	Triton X-100 was added to the conditioned medium to a final concentration of 0.5%. This procedure has been shown to inactivate HIV-1 completely.
Storage:	Store at -70°C
Purity:	More than 95% pure when analyzed by SDS-PAGE

ABL, INC.

9800 Medical Center Drive
Building D
Rockville, Maryland 20850

Phone: 301-816-5200
Fax: 301-816-5438
Email: info@ablinc.com
www.ablinc.com