



HIV-1₄₅₁ gp120

ADVANCED BIOSCIENCE LABORATORIES, INC.



| | |
|-------------------------|---|
| Catalog Number: | 5424 |
| Quantity: | 115 µg |
| Source: | Purified by affinity chromatography and HPLC, from inactivated, concentrated, conditioned medium of HUT78 (clone 6D5) cells chronically-infected with HIV-1 ₄₅₁ isolate. |
| Inactivation Procedure: | Triton X-100 was added to the conditioned medium to a final concentration of 0.5%. This procedure has been shown to inactivate HIV-1 completely. |
| Storage: | Store at -70°C |
| Purity: | More than 95% pure when analyzed by SDS-PAGE |

ADVANCED
BIOSCIENCE
LABORATORIES, INC.

5510 Nicholson Lane
Kensington, Maryland 20895

Phone: 301-816-5200
Fax: 301-816-5438
Email: info@ablinc.com
Website: www.ablinc.com

[Advanced BioScience Laboratories, Inc.](http://www.ablinc.com)