



HIV-1₄₅₁ gp160

ADVANCED BIOSCIENCE LABORATORIES, INC.



Catalog Number:	5427
Quantity:	100 µg
Source:	The conditioned medium of HIV-1 ₄₅₁ infected cells. The gp160 is secreted in the medium. The protein is an oligomer and is probably truncated near the transmembrane domain and runs as a 140 kDa protein in SDS-PAGE.
Inactivation Procedure:	The medium was inactivated with 0.5% Triton X-100 which has been shown to inactivate HIV-1.
Purity:	Purified by immunoaffinity chromatography using mouse monoclonal antibody to HIV-1 gp41. This preparation is free of endotoxin.
Storage:	Store at -70°C

ADVANCED
BIOSCIENCE
LABORATORIES, INC.

5510 Nicholson Lane
Kensington, Maryland 20895

Phone: 301-816-5200
Fax: 301-816-5438
Email: info@ablinc.com
Website: www.ablinc.com

[Advanced BioScience Laboratories, Inc.](http://www.ablinc.com)